

USER GUIDE: Immobilized Metal Ion Affinity Chromatography (IMAC)

IMAC utilizes the surface bonds between the proteins and the histidine residues (or Trp and Cys) binding with divalent metal ions (e.g. Ni^{2+} , Cu^{2+} , Zn^{2+} , Co^{2+}) immobilized via a chelating ligand.

This material is for EBL mam-50 series products induced **His_IDA resin**, **His_NTA resin**, **His_Super resin**, and **IMAC resin**.

Purifying His-tagged recombinant proteins by immobilized metal affinity chromatography (IMAC) continues to grow in popularity. Nickel (Ni^{2+}) is the most commonly used metal ion in IMAC purifications. Both Ni-IDA resin and Ni-NTA resin have higher binding capacities for histidine than other resins with cobalt or other divalent cations. EBL mam-50 series provides different type resin available for the changeable samples.

EBL mam-50 His_Super has lower metal ion leaking, can be compatible with solutions containing EDTA or other chelators and reducers and let us directly purify target protein without dialyzing or buffer changing.

EBL AM-50 series use polyacrylate medium and fully optimize the arm technology for coupling ligands to this matrix. By this, an excellent combination of high binding capacity, high performance and high flow rate addressed. Furthermore, leakage of Ni^{2+} ions is minimized.

Key Features:

- An excellent mechanical and chemical stability.
- Extremely high dynamic binding capacity (DBC), DOUBLED comparing to the products provided by other manufacturer.
- High protein binding capacity and minimal leakage of Ni^{2+} ions
- Extremely stable in high and low pH with excellent mechanical property.
- Compatible with a very wide range of reducing agents, detergents, and other additives
- EBL mam-50 His_Super has high tolerability to EDTA (up to 0.1M) and reduct reagents such as 2-Mer (up to 30mM)

Specifications:

Particle Size (μm)	$50 \pm 5 \%$
Pore Size (\AA)	800
Bead Matrix	Highly cross-linked polyacrylate polymer
Ligand Density	$\approx 150 \mu\text{mol Ni}^{2+}$
Dynamic Binding Capacity (g/ml)	$\approx 90 \text{ mg His-tagged protein/ml medium}$
Max. Linear Flow Rate (cm/h)	≤ 1800
Suggested Flow Rate (cm/h)	150-700
pH Stability	3-13
Pressure (MPa)	10
Chemical Stability	0.01 M HCl, 0.5 M NaOH. Sustained for 1 week at 40 °C. 1 M NaOH, 70 % acetic acid. Sustained for 12 h. 2 % SDS. Sustained for 1 h. 30 % 2-propanol. Sustained for 30 min.
Storage buffer	20 % ethanol
Storage temperature	+4- 28 °C

Quick start protocol:

Preparation of cleared E. coli lysates under native conditions

Notes before starting:

Prepare a bacterial cell pellet (frozen at $\sim 20^\circ\text{C}$ or $\sim 80^\circ\text{C}$)

Lysis may include **Nucleic acid-Cut** (CAT NO: C300102, equivalent to Benzonase nuclease).

Procedure

1. Thaw cell pellet for 15 min on ice. Re-suspend cells in lysis buffer with 10 mM **imidazole** (CAT NO: C300103) at 2~5 ml per gram wet weight.

Note: If the tagged protein does not bind under these conditions, imidazole should be reduced to 1~5 mM. For proteins exhibiting high binding affinities, the imidazole concentration can be increased to 20 mM.

2. Add **lysozyme** (CAT NO: C300106) to 1 mg/ml and **Nucleic acid-Cut** (5 units per ml

of original cell culture volume processed) and incubate on ice for 10 min. Alternative lysis methods can be used, such as ultrasonic, etc.

3. Centrifuge lysate with 10,000 g for 20~30 min at 4 °C to pellet the cellular debris. Supernatant is collected and kept on ice.

4. Any insoluble material must be solubilized before purification under denaturing conditions. Add 5 µl 2 × SDS-PAGE sample buffer to 5µl cleared lysate supernatant for SDS-PAGE analysis.

Purification of 6 × His-tagged proteins under native conditions

Notes before starting

During centrifugation in step 3 of “Preparation of cleared E. coli lysates under native conditions” (see above), equilibrate the Ni-NTA matrix.

Procedure :

1. Pipet 1 ml of Ni-NTA slurry (0.5 ml bed volume) to a 15 ml tube and briefly centrifuge. Remove supernatant and add 2 ml of lysis buffer. Mix by gently inverting. Repeat the centrifugation step and remove the supernatant.

2. Add 4 ml cleared lysate (from step 3 above) to this equilibrated matrix and mix gently by shaking (200 rpm on a rotary shaker) at 4 °C for 60 min.

3. Load the lysate/Ni-NTA mixture into a column with the bottom outlet capped.

4. Remove bottom cap and collect the column flow-through. Save flow through for SDS-PAGE analysis.

5. Wash twice with 5 x bed volume (2.5 ml) wash buffer. Collect wash fractions for SDS-PAGE analysis.

6. Elute the protein 4 times with 0.5 ml elution buffer. Collect the elutions into four tubes and perform SDS-PAGE analysis.

Preparation of cleared E. coli lysates under denaturing conditions

1. Thaw the cell pellet for 15 min on ice and re-suspend in Lysis buffer A or buffer B at 5 ml per gram wet weight.
2. Cells can be lysed in either 6 M Gua-HCl (buffer A) or 8 M Urea (buffer B). It is preferable to lyse the cells in the milder denaturant Urea, so that the cell lysate can be analyzed directly by SDS-PAGE. Gua-HCl is a more efficient solubilization and cell lysis reagent, however, and may be required to solubilize some proteins. Prior to SDS-PAGE analysis, samples containing guanidine must be treated as described in the appendix.
3. The amount of cells required depends on the expression level of the 6 × His-tagged protein and the expression system used. The binding capacity of Ni-resins is protein-dependent and the product provided by EBL normally lies between 40~90 mg/ml medium, doubled comparing to the products provided by other manufacturers.
4. Stir cells for 15~20 min at room temperature or lyse them by gently vortexing, taking care to avoid foaming. Lysis is completed when the solution becomes translucent.
5. Centrifuge lysate with 10,000 × g for 20~30 min at room temperature to pellet the cellular debris. Save supernatant (cleared lysate).
A certain proportion of the cellular protein, including the 6 × His-tagged protein, may remain insoluble and will be located in the pellet. For more complete recovery of the tagged protein, this material must be solubilized as described in this Protocol before purification under denaturing conditions.
6. Add 5 μl 2 × SDS-PAGE sample buffer to 5 μl supernatant and store at ~20°C for future SDS-PAGE analysis.
7. Proceed to protocols for purification under denaturing conditions below.

Purification of 6 x His-tagged proteins from E. coli under denaturing conditions

1. One ml of the EBL Ni-NTA slurry is added to 4 ml lysate and mix gently by shaking (e.g. 200 rpm on a rotary shaker) for 15~60 min at room temperature.

The amount of lysate required depends on the expression level of the 6 × His-tagged protein and the expression system used. The binding capacity of EBL NTA and IDA resins is protein-dependent and normally lies between 40~90 mg/ml medium. For proteins that are expressed at very high levels (50~100 mg of 6 × His-tagged protein per liter of cell culture), a 5 × concentrated cell lysate (resuspend the pellet from a 20 ml culture in 4 ml buffer B) can be used. 4 ml of a 5 × concentrated cell lysate in buffer B will contain approximately 1~2 mg of 6 × His-tagged protein. For much lower expression levels (1~5 mg/liter), 200 ml of cell culture should be used for a 50 × concentrated cell lysate (4 ml cell lysate = 0.2~1 mg of 6 × His-tagged protein).

2. Load lysate–resin mixture carefully into an empty column with the bottom cap still attached.

3. Remove the bottom cap and collect the flow-through. Collect flow-through for SDS-PAGE analysis.

4. Wash twice with 4 ml buffer C. Keep washing fractions for SDS-PAGE analysis.

5. Elute the recombinant protein 4 times with 0.5 ml buffer D, followed by 4 times with 0.5 ml buffer E. Collect fractions and analyze by SDS-PAGE. Monomers generally elute in buffer D, while multimers, aggregates, and proteins with two 6 × His tags will generally elute in buffer E.

General suggestions

- ※ Keep cells and protein solutions at 0–4 °C at all times to prevent protein degradation. The addition of protease inhibitors may also be necessary.
- ※ Low concentrations of imidazole in lysis and wash buffers minimize nonspecific binding and reduce the amount of contaminating proteins.
- ※ The capacity of this product is doubled comparing with those provided by other suppliers. Because of this, more imidazole should be needed to complete elution of the target protein. For this purpose, you can use the following strategies or combination of them: **I.** more volumes of the elution buffer, or **II.** to increase the concentration of imidazole in elution buffer. **III.** An adjustment of elution pH buffer can be very effective in many cases. The pH should be optimized via assays

- ※ Elution conditions are highly reproducible, but optimum conditions must be determined for each His-Tag fusion protein purified. Although several elution methods are effective (imidazole, pH, and EDTA), imidazole is the mildest and is recommended under native conditions, when the protein would be damaged by a reduction in pH, or when the presence of metal ions in the elution may have an adverse effect on the purified protein. Monomers generally elute at approximately pH 5.9, whereas aggregates and proteins that contain more than one His-Tag sequence elute at approximately pH 4.5.

Appendix: Buffers recipe

Buffers for purification under **native conditions** (from *E. coli* and insect cells)

--Lysis buffer (1 liter):

50 mM NaH₂PO₄, 300 mM NaCl, 10 mM imidazole, adjust pH to 8.0 using NaOH.

--Wash buffer (1 liter):

50 mM NaH₂PO₄, 300 mM NaCl, 20 mM imidazole, adjust pH to 8.0 using NaOH.

--Elution buffer (1 liter):

50 mM NaH₂PO₄, 300 mM NaCl, 250 mM imidazole, adjust pH to 8.0 using NaOH.

Buffers for purification under **denaturing conditions**

--Lysis buffers A (1 liter):

100 mM NaH₂PO₄, 10 mM Tris, 6 M Gua-HCl, adjust pH to 8.0 using NaOH.

--Lysis buffer B (1 liter):

100 mM NaH₂PO₄, 10 mM Tris, 8 M urea, adjust pH to 8.0 using NaOH.

--Buffer C (1 liter):

100 mM NaH₂PO₄, 10 mM Tris, 8 M Urea, adjust pH to 6.3 using HCl.

--Buffer D (1 liter):

100 mM NaH₂PO₄, 10 mM Tris, 8 M Urea, adjust pH to 5.9 using HCl.

--Buffer E (1 liter):

100 mM NaH₂PO₄, 10 mM Tris, 8 M Urea, adjust pH to 4.5 using HCl.

Note: Due to the dissociation of urea, the pH of Buffers B, C, D, and E should be adjusted immediately prior to use. Do not autoclave.

Buffers for purification from mammalian cells under native conditions

PBS: 50 mM potassium phosphate, pH 7.2; 150 mM NaCl

Lysis buffer (1 liter):

50 mM NaH₂PO₄, 300 mM NaCl, 10 mM Imidazole, 0.05 % Tween 20, adjust pH to 8.0 using NaOH.

Wash buffer (1 liter):

50 mM NaH₂PO₄, 300 mM NaCl, 20 mM imidazole, 0.05 % Tween 20, adjust pH to 8.0 using NaOH.

Elution buffer (1 liter):

50 mM NaH₂PO₄, 300 mM NaCl, 250 mM imidazole, 0.05 % Tween 20, adjust pH to 8.0 using NaOH

Preparation of guanidine-containing samples for SDS-PAGE

Since fractions that contain Gua-HCl will form a precipitate when treated with SDS, they must either be diluted with H₂O (1:6), dialyzed before analysis, or separated from the Gua-HCl by TCA precipitation.

TCA precipitation:

1. Dilute 10–25 µl protein samples to 100 µl with ddH₂O
2. Add 100 µl of 10 % TCA.
3. Leave on ice for 20 min; centrifuge for 15 min in a microcentrifuge.
4. Wash pellet with 100 µl of ice-cold ethanol, dry, and resuspend in sample buffer.

If there is any trace of Gua-HCl present, samples must be loaded immediately after heating for 7 min at 95 °C.

Reuse of Ni-NTA Resin:

The reuse of EBL mam-50 series resin depends on the nature of the sample and should only be performed with identical recombinant proteins. We recommend a maximum of 5 runs.

If the resin changes from light blue to brownish-gray, the following regeneration procedure is recommended.

Wash the resin / column successively with:

1. 2 vol regeneration buffer (6 M Gua-HCl, 0.2 M acetic acid).
2. 5 vol water.
3. 3 vol 2% SDS.
4. 1 vol 25 % ethanol.
5. 1 vol 50 % ethanol.
6. 1 vol 75 % ethanol.
7. 5 vol 100 % ethanol.
8. 1 vol 75 % ethanol.
9. 1 vol 50 % ethanol.
10. 1 vol 25 % ethanol.
11. 1 vol water.
12. 5 vol 100 mM EDTA, pH 8.0.
13. 10 vol water.
14. Recharge the column with 2 vol 100 mM NiSO₄.
15. Wash the column with 2 vol water.
16. Wash the column with 2 vol regeneration buffer.
17. Equilibrate with 2 vol of a suitable buffer (e.g. Lysis Buffer).

Ordering information:

CAT NO:	Product	Quantity
pct-m500115	His_IDA resin prepacking colume	1mlx1/ 1mlx5/ 5mlx1
mam-5001005	MAM-50 His_IDA resin	5ml/ 20ml/ 100ml
pct-m500215	IMACA resin prepacking colume	1mlx1/ 1mlx5/ 5mlx1
mam-5002005	MAM-50 IMAC resin	5ml/ 20ml/ 100ml
pct-m500315	His_NTA resin prepacking colume	1mlx1/ 1mlx5/ 5mlx1
mam-5003005	MAM-50 His_NTA resin	5ml/ 20ml/ 100ml
pct-m500515	His_Co++ resin prepacking colume	1mlx1/ 1mlx5/ 5mlx1
mam-5005005	MAM-50 His_Co++ resin	5ml/ 20ml/ 100ml
pct-m500615	His_Super resin prepacking colume	1mlx1/ 1mlx5/ 5mlx1
mam-5006005	MAM-50 His_Super resin	5ml/ 20ml/ 100ml

Related product:

CAT NO:	Product	Quantity
C300102	Nucleic acid-Cut	5KU/ 25KU
C300103	imidazole	100g/ 500g
C300106	lysozyme	10g/ 50g
C103500	Tris base	500g